

# Configuration Planning

## ENDOMAT® SELECT

The Choice is Yours



# ENDOMAT® SELECT

The Choice is Yours

ENDOMAT® SELECT is a single roller pump for the irrigation or suction of fluids during surgical and diagnostic procedures.

The pump can be equipped with various software modules so that the unit can be used for multiple disciplines in various rooms according to the respective indication.

The unit offers various software packages which the user can select and combine according to individual needs.




- UP210 **ENDOMAT® SELECT SCB**, with integrated SCB module, suction and irrigation pump, incl. mains cord, power supply 100-240 VAC, 50/60 Hz

When placing an initial order, please tick the box for UP210. Furthermore, at least one software module must be included in the order. Please tick the appropriate box in the document.

Please do not tick the box if a software upgrade for an existing unit is requested.

The following software packages are available for ENDOMAT® SELECT:

Clicking on one of the fields highlighted in color takes you to the relevant page.



The following tubing sets are available for ENDOMAT® SELECT:



For pressure-controlled applications such as, for example, hysteroscopy, cystoscopy or interventions in arthroscopy or spine surgery:

031523-10\*


**Tubing Set**, irrigation, pressure-controlled, STERILE , sterile, for single use, package of 10, for use with KARL STORZ HAMOU® ENDOMAT® 26331120-1 and ENDOMAT® SELECT UP210

Procedures that use a pressure-controlled tubing set bear the adjacent mark.

## SURGERY

For flow-controlled applications such as, for example, laparoscopy, thoracoscopy, proctology:




031524-10\* **Tubing Set**, irrigation, flow-controlled, STERILE , sterile, for single use, package of 10, for use with KARL STORZ HAMOU® ENDOMAT® 26331120-1 and ENDOMAT® SELECT UP210



Procedures that use a flow-controlled tubing set bear the adjacent mark.

For hysteroscopic tissue ablation with the IBS® shaver system or the removal of tissue chips during resection (TUR-B/P):




030647-10\* **Tubing Set**, suction, STERILE , direct suction, sterile, for single use, package of 10, for use with KARL STORZ ENDOMAT® SELECT UP210



Procedures that use a tubing set for direct suction bear the adjacent mark.

For use in combination with the CALCUSON lithotripsy unit in urology:



031647-10\* **Tubing Set**, suction, STERILE , via bottle, sterile, for single use, package of 10, for use with KARL STORZ ENDOMAT® SELECT UP210



Procedures requiring a suction bottle bear the adjacent mark.

For use with the CLEARVISION® functionality in ENT/NEURO:



031529-10\* **Tubing Set**, irrigation, CV, sterile, for single use, package of 10, for use with KARL STORZ ENDOMAT® SELECT UP210




Procedures that use a CLEARVISION® functionality bear the adjacent mark.


In the daily surgical routine in the OR, the ENDOMAT® SELECT delivers a powerful yet low-noise irrigation of fluids during interventions related to the following specialties:


- Laparoscopy
- Thoracoscopy
- Proctology

### Software:

- UP601 **SURGERY Software**, license, allows selection of the procedures “LAP”, “THOR” and “PROCTO”
- Promotion Pack (031524-03\*)** STERILE , **Tubing Set**, irrigation, FC, sterile, for single use, package of 3

### Further accessories:

- 031524-10\* **Tubing Set**, irrigation, flow-controlled, STERILE , sterile, for single use, package of 10, for use with KARL STORZ ENDOMAT® HAMOU® 26331120-1 and ENDOMAT® SELECT UP210



## (OFFICE) HYSTEROSCOPY

ENDOMAT® SELECT provides pressure-controlled dilation of the uterine cavity for diagnostic and operative hysteroscopy in the office setting or the clinical outpatient department.

### Software:

<input type="checkbox"/>	UP602	<b>HYSTEROSCOPY Software</b> , license, allows selection of the procedure "HYS"	
<input checked="" type="checkbox"/>		<b>Promotion Pack (031523-03*)</b>	<b>Tubing Set</b> , irrigation, PC, sterile, for single use, package of 3

### Further accessories:

031523-10\* **Tubing Set**, irrigation, pressure-regulated, sterile, for single use, package of 10, for use with KARL STORZ ENDOMAT® HAMOU® 26331120-1 and ENDOMAT® SELECT UP210

## Hysteroscopic Tissue Ablation with the Intrauterine BIGATTI Shaver IBS®

In conjunction with the IBS® shaver, the ENDOMAT® SELECT can be used for the controlled suction of resected polyps, myomas or placental remnants.



### Software:

<input type="checkbox"/>	UP603	<b>IBS® SHAVER</b> , software, license, allows selection of the procedure "IBS®"	
<input checked="" type="checkbox"/>		<b>Promotion Pack (030647-03*)</b>	<b>Tubing Set</b> , suction, DS, sterile, for single use, package of 3

### Mandatory accessories:

030647-10\* **Tubing Set**, suction, direct suction, sterile, for single use, package of 10, for use with KARL STORZ ENDOMAT® SELECT UP210

20701070 **Control Cable**, connectors 2x LEMO 5-pol 0°, length 100 cm, for transmission of footswitch control signal between ENDOMAT® SELECT UP210 and UNIDRIVE® S III 20701020



## UROLOGY

ENDOMAT® SELECT accommodates the wide range of applications in the field of urology, allowing seamless integration into existing systems and providing pressure-regulated irrigation for diagnostic and operative interventions.

### Software:

<input type="checkbox"/>	UP604	<b>UROLOGY Software</b> , license, allows selection of the procedures "CYST", "RES", "URS", "CALCUSON", and "PCN"	
<b>Please select the required Promotion Pack:</b>			<b>STERILE</b>
<input type="checkbox"/>		<b>Promotion Pack (031523-03*) Tubing Set</b> , irrigation, PC, sterile, for single use, package of 3	
or	<input type="checkbox"/>	<b>Promotion Pack (030647-03*) Tubing Set</b> , suction, DS, sterile, for single use, package of 3	
or	<input type="checkbox"/>	<b>Promotion Pack (031647-03*) Tubing Set</b> , suction, BS, sterile, for single use, package of 3	

### Further accessories:


031523-10\* **Tubing Set**, irrigation, pressure-controlled, sterile, for single use, package of 10, for use with KARL STORZ HAMOU® ENDOMAT® 26331120-1 and ENDOMAT® SELECT UP210

030647-10\* **Tubing Set**, suction, direct suction, sterile, for single use, package of 10, for use with KARL STORZ ENDOMAT® SELECT UP210

031647-10\* **Tubing Set**, suction, via bottle, sterile, for single use, package of 10, for use with KARL STORZ ENDOMAT® SELECT UP210

20014130 **One-Pedal Footswitch**, digital, one-stage


20701070 **Control Cable**, connectors 2x LEMO 5-pol 0°, length 100 cm, for transmission of footswitch control signal between ENDOMAT® SELECT UP210 and UNIDRIVE® S III 20701020 such as between ENDOMAT® SELECT UP210 and CALCUSON 27610020




## ARTHROSCOPY


In the field of arthroscopy, the ENDOMAT® SELECT regulates pressure-controlled irrigation of fluids to maintain joint distension.

### Software:

- UP605 **ARTHROSCOPY Software**, license, allows selection of the procedures “KNEE”, “HIP”, “SHOULDER” and “SMALL JOINTS”,  
for use with ENDOMAT® SELECT UP210
- Promotion Pack (031523-03\*)** STERILE   
**Tubing Set**, irrigation, PC, sterile, for single use, package of 3

### Further accessories:


- 031523-10\* **Tubing Set**, STERILE   
irrigation, pressure-regulated, sterile, for single use, package of 10, for use with KARL STORZ HAMOU® ENDOMAT® 26331120-1 and ENDOMAT® SELECT UP210
- 20014130 **One-Pedal Footswitch**, digital, one-stage
- UP006 **Control Cable**,  
connectors 2x LEMO 5-pol 0°,  
length 100 cm, for transmission of footswitch control signal between ENDOMAT® SELECT UP210 and UNIDRIVE® S III ARTHRO SCB 28723020-1




## ENT/NEURO


ENDOMAT® SELECT includes the following software module for irrigation of the lens

### Software:

- UP606 **ENT/NEURO Software**, allows selection of the procedure “CLEARVISION®”
- Promotion Pack (031529-03\*)** STERILE   
**Tubing Set**, irrigation, CV, sterile, for single use, package of 3, for use with KARL STORZ ENDOMAT® SELECT UP210

### Mandatory accessories:

- 031529-10\* **Tubing Set**, STERILE   
irrigation, CV, sterile, for single use, package of 10, for use with KARL STORZ ENDOMAT® SELECT UP210
- 20014230 **One-Pedal Footswitch**, digital, two-stages




## SPINE SURGERY

ENDOMAT® SELECT offers the following software module for the irrigation of fluids during interventions in spine surgery:

### Software:

- |                                     |                                    |   |
|-------------------------------------|------------------------------------|---|
| <input type="checkbox"/>            | UP607                              | <b>SPINE Software</b> , license, allows selection of the procedures "LUMBAR" and "THORACIC" |
| <input checked="" type="checkbox"/> | <b>Promotion Pack (031523-03*)</b> | <b>Tubing Set</b> , irrigation, PC, sterile, for single use, package of 3                   |



### Further accessories:

- |            |  |
|------------|--|
| 031523-10* | <b>Tubing Set</b> ,<br>irrigation, pressure-controlled, sterile, for single use,<br>package of 10, for use with KARL STORZ HAMOU®<br>ENDOMAT® 26331120-1 and ENDOMAT® SELECT UP210 |
|------------|--|

## ADVANCED

In order to activate extended settings and additional functions, the "Advanced Package" can be selected as an additional\* software package.

### Software:

- |                          |       |   |
|--------------------------|-------|---|
| <input type="checkbox"/> | UP610 | <b>Advanced Package</b> , software, license |
|--------------------------|-------|---|

\* The advanced package only functions when at least one other software package is available.

The following features are available with the installation of ADVANCED Software UP610:

- Adjustment of pressure (100/300/500 mmHg) during SURG applications
- Adjustment of inflow up to 3500 ml/min during SURG applications
- Adjustment of inflow (200/400/600 ml/min) during HYS, SPINE & URO applications
- Adjustment of inflow (1500/2000/2500 ml/min) during ARTHRO applications
- Duration of BOOST during ARTHRO & URO applications
- Adjustment of pump position (height) to improve pump performance
- Adjustment of volume of the button tones
- Adjustment of brightness
- SCB Integration

The logo features the number '75' in a large, bold, blue sans-serif font. The '5' has a white circle in its center, resembling an endoscopic lens. Below the '75', the word 'Years' is written in a blue, cursive script font.

*Shaping the Future  
of Endoscopy with you*

**STORZ**  
**KARL STORZ—ENDOSKOPE**

*THE DIAMOND STANDARD*

KARL STORZ SE & Co. KG  
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany  
Postfach 230, 78503 Tuttlingen/Germany  
Telefon: +49 7461 708-0  
Telefax: +49 7461 708-105  
E-Mail: [info@karlstorz.com](mailto:info@karlstorz.com)  
[www.karlstorz.com](http://www.karlstorz.com)

