

Enseal®

ENSEAL® G2 Curved and Straight Tissue Sealers

Secure sealing and versatility
in the palm of your hand.



ETHICON
PART OF THE *Johnson & Johnson* FAMILY OF COMPANIES

Shaping
the future
of surgery

Secure **sealing** and **versatility** in the palm of your hand



Secure

- Twice the compression at the distal tip of the jaw than LigaSure™ 5mm¹
- Less proximal end bleeding† and more consistent seal strength over multiple activations‡ than LigaSure™ 5mm
- Seals vessels up to 7mm diameter and lymphatics


Versatile

- Delivers energy to tissue even with the jaws open, enabling spot coagulation, enhancing dissection and creating otomies


† In a preclinical study utilizing an animate porcine mesentery model, ENSEAL® G2 Curved and Straight had significantly less bleeding at the proximal end than LigaSure™ 5mm (0.0% vs. 24.0%), P=0.001. 062106-220324

‡ In a benchtop study with ENSEAL® G2 Curved and Straight using 5mm porcine carotid arteries, vessel burst pressure was measured over 52 activation cycles P<0.05. 094215-210708

¹ In side-by-side benchtop testing, average compression of jaws used in ENSEAL® G2 Straight was 169 psi (standard deviation 61 psi) vs. 84 psi (standard deviation 4 psi) for LigaSure™ 5mm, p<0.0001. 093419-210626



ENSEAL® G2 Tissue Sealers vs LigaSure™



Secure **sealing** and **versatility**

A comprehensive redesign – in both the handle and activation sequence – ensures comfortable and confident operation.

Integrated activation


- Provides energy activation prior to cutting

Audible tones

- Provides feedback throughout the seal cycle

Ergonomic design

- Designed for comfort and performance



The Ethicon Endo-Surgery Generator
Revolutionary power
for the ENSEAL® G2 Tissue Sealers (sold separately)



ENSEAL® G2 Tissue Sealers

The complete ENSEAL® G2 Tissue Sealer Portfolio

All ENSEAL® and HARMONIC® devices are compatible with the Ethicon Endo-Surgery Generator.



ENSEAL® G2 Curved

Code	Shaft diameter	Shaft length	Qty
NSLG2C14	5mm	14cm	6
NSLG2C25	5mm	25cm	6
NSLG2C35	5mm	35cm	6
NSLG2C45	5mm	45cm	6



ENSEAL® G2 Straight

Code	Shaft diameter	Shaft length	Qty
NSLG2S14	5mm	14cm	6
NSLG2S25	5mm	25cm	6
NSLG2S35	5mm	35cm	6
NSLG2S45	5mm	45cm	6

ENSEAL® G2 Curved and Straight Tissue Sealers – Secure sealing and versatility to do more.

To experience the ENSEAL® G2 Tissue Sealers, contact your Ethicon Sales Professional
Call 1-877-ETHICON. Connect at www.ethicon.com/ensealg2

For complete indications, contraindications, warnings, precautions, and adverse reactions, please reference full package insert.

ETHICON | Shaping
PART OF THE *Johnson & Johnson* FAMILY OF COMPANIES | the future
of surgery

Enseal®