

Designed to deliver reliable performance on an easy to use platform



ENDOPATH™ ETS-Flex 45

Compression before firing

- Designed to pre-compress tissue for consistent staple formation

Ease of use

- Designed to clamp, staple, cut and release all with one hand
- Designed with 5 angle increments of articulation to each side enabling a more precise approach

ETHICON


Johnson & Johnson SURGICAL TECHNOLOGIES

Product components



ENDOPATH™ ETS Articulating Linear Cutter 45mm*

CODE	LENGTH	ARTICULATION	QUANTITY PER SALES UNIT
ATS45	340mm	Yes	3
ATS45NK	340mm	Yes	3



ENDOPATH™ ETS Reloads

CODE	COLOR	STAPLE ROWS	TISSUE THICKNESS	OPEN STAPLE HEIGHT	CLOSED STAPLE HEIGHT	QUANTITY PER SALES UNIT
TR45W	White	6	Vascular/Thin	2.5mm	1.0mm	12
6R45B	Blue	6	Regular	3.5mm	1.5mm	12
TR45G	Green	6	Thick	4.1mm	2.0mm	12

How to order

Electronic ordering options

Purchases can be made through your distributor or direct at Johnson & Johnson Health Care Systems, Inc. (JJHCS).

The following electronic order placement methods are preferred:

- J&J Gateway—www.jnjgateway.com (for questions about your order, visit the Web site or call 1-866-JNJ-GATE)
- Electronic Data Interchange—JJHCS Help Line: 1-800-262-2888

Nonelectronic/manual ordering options

JJHCS—Call 1-800-255-2500 between the hours of 8:30 am and 8:00 pm Eastern time, or fax your order to 1-732-562-2212.

Customer support

For product use assistance, clinical guidelines, service and repair, emergency assistance, copy of a 510(k) clearance letter, or complaints, please contact our Customer Support Center by calling 877-ETHICON (384-4266). Our support center is staffed 24 hours a day, 7 days a week by qualified nurses to answer your product-related questions.

Visit www.ethicon.com for more information about the ENDOPATH ETS Articulating Linear Cutter.

For complete indications, contraindications, warnings, precautions, and adverse reactions, please reference full package insert.

ETHICON

Johnson & Johnson SURGICAL TECHNOLOGIES