

Single-use Biopsy Forceps

EndoJaw

Comprehensive Lineup Delivering Optimal Performance.



Deliver Optimal Performance

Olympus offers an extensive range of single-use biopsy forceps. Featuring superior construction and delivering reliable performance, EndoJaw devices are designed to work seamlessly with a broad array of Olympus endoscopes.

Comprehensive Lineup

EndoJaw biopsy forceps are available in a wide variety of configurations, including alligator jaw and needle designs as well as a choice of diameters: 2.0 mm, 2.8 mm, 3.2 mm, and 3.7 mm.

Alligator jaw with needle
3.7 mm

Alligator jaw
3.2 mm

EndoJaw Large
Alligator jaw with needle
2.8 mm

Oval jaw with needle
2.8 mm

Oval jaw
2.0 mm

Simplified Targeting

Proprietary swing jaw design mechanism facilitates the positioning of the forceps cups at biopsy sites and allows for tangential approaches.



Pristine Precision


The EndoJaw's cups are made of stainless steel for optimal precision.

Smooth Insertion Sheath

All EndoJaw models are equipped with a sheath design that provides optimal flexibility to maximize insertion capability into the endoscope channel.

Double-step Blade Performance

The unique dual-thickness design of the EndoJaw's cups provides clean specimen retrieval by use of sharp cutting performance with thinner edges, while the thicker body ensures durability and reliability.



Position forceps.

Grasp tissue.

Cut biopsy.

EndoJaw

Single-use Biopsy Forceps

| Specifications for Upper GI | | | | | | |
|-----------------------------|----------------------|---------------|-------------|----------|----------------|--------------|
| Catalog number | Minimum channel size | Cup type | Feature | Quantity | Working length | Sheath color |
| FB-210K.A | 2.8 mm | Alligator Jaw | - | 20/box | 1550 mm | Gray |
| FB-220K.A | | | With needle | | | |
| FB-230K.A | | Oval Jaw | - | | | |
| FB-240K.A | | | With needle | | | |
| FB-211K.A | 2.0 mm | Alligator Jaw | - | | | |
| FB-221K.A | | | With needle | | | |
| FB-231K.A | | Oval Jaw | - | | | |
| FB-241K.A | | | With needle | | | |

| Specifications for Lower GI | | | | | | |
|-----------------------------|-----------------------------|---------------|-------------|----------|----------------|--------------|
| Catalog number | Minimum channel size | Cup type | Feature | Quantity | Working length | Sheath color |
| FB-210U.A | 2.8 mm Standard Capacity | Alligator Jaw | - | 20/box | 2300 mm | Blue |
| FB-220U.A | | | With needle | | | |
| FB-230U.A | | Oval Jaw | - | | | |
| FB-240U.A | | | With needle | | | |
| FB-215U | 2.8 mm Large Capacity | Alligator Jaw | - | | | |
| FB-225U | | | With needle | | | |
| FB-235U | | Oval Jaw | - | | | |
| FB-245U | | | With needle | | | |
| FB-214U | 3.2 mm | Alligator Jaw | - | 20/box | | |
| FB-224U | | | With needle | | | |
| FB-234U | | Oval Jaw | - | | | |
| FB-244U | With needle | | | | | |
| FB-212U | 3.7 mm | Alligator Jaw | - | | 5/box | |
| FB-222U | | | With needle | | | |
| FB-232U | | Oval Jaw | - | | | |
| FB-242U | | | With needle | | | |

| Specifications for Pulmonary | | | | | | |
|------------------------------|----------------------|---------------|-------------|----------|----------------|--------------|
| Catalog number | Minimum channel size | Cup type | Feature | Quantity | Working length | Sheath color |
| FB-211D.A | 2.0 mm | Alligator Jaw | - | 20/box | 1150 mm | Green |
| FB-221D.A | | | With needle | | | |
| FB-231D.A | | Oval Jaw | - | | | |
| FB-241D.A | | | With needle | | | |

Specifications, design, and accessories are subject to change without any notice or obligation on the part of the manufacturer.



OLYMPUS EUROPA SE & CO. KG
 Postbox 10 49 08, 20034 Hamburg, Germany
 Wendenstrasse 14-18, 20097 Hamburg, Germany
 Phone: +49 40 23773-0, Fax: +49 40 233765
www.olympus-europa.com